


ZIC®-HILIC Separation of Sugar Phosphates

Chromatographic Conditions

Column: ZIC®-HILIC, PEEK 150 x 4.6 mm, 5 µm, 200 Å (P/N 150455)
 Injection: 1 µL injection of a standard solution diluted in in acetonitrile/water (60:40)
 1: D-Ribulose 1,5-bisphosphate sodium salt 0.07%(w/v)
 2: 2,3-Diphospho-D-glyceric acid penta sodium salt 0.15%(w/v)
 Detection: UV @ 210 nm (UFS 1.0 V)
 Column Temperature: 40 °C
 Flow Rate: 1.0 mL/min
 Mobile Phase (v/v): Acetonitrile/ 40 mM Potassium phosphate pH 3.0 (65:35 v/v)
 (total ionic strength 14 mM)



Chromatographic Data

No.	Compound	Time (min)	Retention factor
	t_0 void volume	1.5	-
1	D-Ribulose 1,5-bisphosphate	6.7	3.5
2	2,3-Diphospho-D-glyceric acid	8.2	4.5

Data kindly provided by Dr Norikazu Nagae, ChromaNik Technologies Inc., Japan